
美容家電のパイオニア
YA-MAN

Instruction Manual

家庭用美容器
Uruoi HOT&COOL

Platinum White
Uruoi HOT&COOL

HB-9-HK



- Before using this product, please read these instructions carefully.
- Refer to the Safety Precautions section for important information.
- Keep this Instruction Manual for future reference.

Safety Precautions

- **Please read all instructions carefully for safe and proper use.**
- **Use only as directed.**
- **Immediately discontinue use if any malfunction or unusual reaction to treatment occurs.**
- **Precautions listed in the instruction manual are very important. For your safety, please follow them carefully.**
- **Precautions are classified as DANGER, WARNING, or CAUTION in order to clearly specify the level of danger and possible injury or damage that may occur if the device is misused.**

DANGER	**Improper use of this device poses a severe and likely risk of death or serious injury.**
WARNING	**Improper use of this device may pose a risk of death or serious injury.**
CAUTION	**Improper use of this device may result in fire, injury, or property damage.**

This symbol shows you what you must NOT do.

This symbol shows you what you must do.

This symbol shows you what you must be careful with.

DANGER

Do not use under any of the following conditions as this may cause malfunction of medical devices.

Do not use in conjunction with any medical electronic equipment such as a pacemaker, oxygenator, electrocardiograph, etc.

WARNING

Improper use of this device may cause injury to your skin or body

- *While pregnant or nursing or during menstruation*
- *While growing or during adolescence*
- *If you suffer from atopic dermatitis or have unusually sensitive skin*
- *If you have skin allergies or bleeding diathesis*
- *If you have skin irritation or redness caused by sunburn*
- *There is suppuration or inflammation of the skin*
- *You experience skin irritation caused by cosmetics*
- *Have itchy pigmentation or pathologic pigmentation*
- *Have viral warts or melanoderma*
- *Have high blood pressure*
- *Have weak skin sensibility*
- *If you are sick or feeling unwell*
- *If you are allergic to metals*
- *If you have capillary dilatation due to long term use of steroid hormonal drugs or a disorder of liver function*

Do not use if you suffer from any of the following conditions, or on body areas listed below as this may cause injury:

- *Acute disorder, febrile disorder, contagion, fissures, feeling unwell, extremely exhausted, tuberculosis, abnormal blood pressure, malignancy, heart disease, or if you have consumed alcohol or medication that causes drowsiness.*
- *Do not use inside the mouth or around the pubic area, mucous membranes, on eyelids or eyeballs, or around the eyes. Do not use over broken bones or areas where you have had plastic surgery.*
- *If you are not sure, consult your doctor before use.*

Use at AC100-240V only to avoid a potential fire or electrical shock hazard.

Do not disassemble, attempt to repair, or modify this device in any way as this may cause a fire, electrical shock hazard, or injury. Consult your local agent if repair is needed.
*Warranty applies only to use within country of purchase. Warranty repair for foreign customers is not available.

Do not use if the AC adapter is damaged or the wall socket is loose.
This may cause electrical shock hazard, short circuit, or fire.

Hold the head of the plug instead of the cord when unplugging the adapter from the wall socket.
Damage to the cord may cause electrical shock, short circuit, or fire hazard.

WARNING

Unplug the AC adapter from the wall socket when not in use or when cleaning as this may cause an electrical shock hazard, short circuit, or fire.

If any problem occurs during use, turn off the power and unplug the AC adapter immediately.
Discontinue use if there is smoke or an unusual noise or odor.

When using the AC adapter, plug it directly into the wall socket.
Do not use an extension cord as this may cause an electrical shock hazard, short circuit, or fire.

Do not use or store near fire, under direct sunlight, or in humid places such as the bathroom.
Do not place near heating equipment as this may cause a fire or explosion hazard.

Do not use on children or allow children to use the device. Keep out of reach of children.
Do not allow disabled persons to use the device unattended.
Avoid use on persons with sensitive skin or on those who cannot express their intentions as this may cause injury.

Do not allow the main unit, AC adapter, or cord to get wet.
Avoid use near water as this may cause an electrical shock hazard, short circuit, or fire.

Do not use near flammable substances such as thinner, benzene, nail polish remover, or hair spray as this may cause fire, electrical shock hazard, or injury.

Do not drop or hit the main unit. Dropping the device may cause malfunction or injury.



WARNING 注意

To avoid accidents, injury, or malfunction:
- Do not use in the bathroom or wet areas
- Do no immerse the device in water or wash it.

CAUTION 注意

Follow the precautions below to avoid accidents, injury, or malfunction:
- Do not apply over the same area for more than 5 seconds.
- Stop use immediately and consult a doctor if there is any unusual skin reaction during or after use.
- Apply according to your particular skin type and condition.
- Avoid use if you are sunburned, or in conjunction with peeling or face mask treatments.
- Always turn off the power after use.

Stop use if treatment is unpleasant or you feel uncomfortable.

To avoid accidents, injury, or malfunction:
- Do not use if you have any unusual skin condition such as rash, irritation, or itching.
- Do not use for any purpose other than as directed.
- Do not use if the device is damaged, has been dropped, or subjected to impact.
- Do not use for professional treatments and avoid continuous or excessive use.

Do not use longer than 10 minutes per treatment:
When used for 10 minutes, the Cool plate may not fully cool. Allow at least a 1 hour interval before next use.

Do not block the vents during use. Do not allow water or any foreign substance inside the vents.

Do not ap p l y to p ets o r animals as this may cause injury.

Use in Foreign Countries

WARNING

Use at AC100-240V only to avoid a potential fire or electrical shock hazard.
*** Warranty only applies to use within the country of purchase. Warranty repair for customers abroad is not available.**

Be sure to check the plug shape and the voltage before using. Using this device with an improper plug or at the incorrect voltage may cause overheating, fire, or malfunction.


NOTE: This product can be operated with 100-240V. If the plug of the supplied AC adapter does not match the plug of your location, use a plug adapter (not provided).


Main plug configurations and countries of use


Type	Plug Configuration	Countries of Use
A		US, Canada, Hawaii, Spain, Guam, Saipan, Thailand, Korea, Taiwan, China, Phillipines, Vietnam, Brasil, Turkey, etc.
B3		UK, Switzerland, China, Hong Kong, India, Singapore, Phillipines, Egypt, South Africa, etc.
BF		Canada, UK, Switzerland, Thailand, Korea, Singapore, China, Hong Kong, Egypt, South Africa, etc.
C		Italy, Holland, Switzerland, Spain, Germany, France, India, Thailand, Korea, China, Hong Kong, Malaysia, Turkey, etc.
O		UK, France, Australia, New Zealand, Korea, Taiwan, China, Phillipines, Argentina, etc.

- The plug configuration and voltage varies depending on the region and the type of accommodation. Always be sure to check for compatibility before using the device.
- For information about plug adapters, consult a travel agency or a consumer electronics retailer.


● Main Unit







● AC Adapter




** Design and specifications of the main unit and accessories are subject to change for improvements without prior notice.*

Indicator Panel

• LED Indicators

Mode and Level are shown on the Indicator Panel (below).



• HOT Lamp

When you press and hold the POWER/MODE button, both the Warm and Cool plate turn on and the HOT lamp flashes red.
When the temperature of the Warm plate reaches 42°C, the HOT lamp will turn solid red.

• MODE Indicator

1 CLEAN mode (Power on) → 2 MOISTURE mode → 3 SMOOTH mode → 4 EYE CARE mode → 5 HOT & COOL mode

When you press and hold the POWER/MODE button, power turns on at 1-CLEAN mode. Press the POWER/MODE button while the power is on to choose the mode. Each time POWER/MODE is pressed, the mode changes as shown at left.

*** 1 CLEAN mode** : Skin cleansing with ion

*** 2 MOISTURE mode** : Moisturizes deep under the skin* with iontophoresis

*** 3 SMOOTH mode** : Moisturizes the skin surface with iontophoresis


*** 4 EYE CARE mode** : Moisturizes around the eyes with iontophoresis

*** 5 HOT/COOL mode** : Applies Warm and Cool treatment

** To stratum corneum.*


Using Modes 1 - 4

Power will automatically turn off after 5 minutes from when the HOT lamp turns solid red. During modes 1-4, a buzzing sound will indicate when 3 minutes has elapsed.




Using Mode 5 : HOT & COOL

Power will automatically turn off after 12 minutes from when the HOT lamp turns solid red.



You may hear a buzzing sound during use of the HOT & COOL mode. This is the cooling system operating and is not a malfunction.

- **Level Indicator**



When you press the LEVEL button while the power is on, the intensity of ion treatment can be selected.

Each time the LEVEL button is pressed, the level will change as shown at left with the solid red lamps indicating the selected level. Choose a level that you find comfortable.


*When beginning treatment, always start at level 1.

*For modes 1-4, the level is adjustable and will always be indicated by solid red lamps.

- **Charging Indicator**

The Charge Lamp will flash red while charging. When charging has completed, the lamp will turn solid red. (Refer to “Charging the battery” on p.11).

* Charge lamp will not turn on during use.



Precautions for Use

Before use, please read this “Precautions for Use” section carefully and operate the device only as instructed.

- **Continuous treatment time**

  Do not use longer than 10 minutes per treatment.
  When repeating a treatment, wait at least 1 hour.
  *Temperature of the plates may become unstable when over-used.

- **Vent**

  Do not block or cover the cooling vents while holding the device.
  Be careful to prevent water or any foreign substance from entering the vent holes.

- **Hand electrode**

  Always hold the device so that your hand is in contact with the hand electrode.
  Ion treatment will not work if you are not touching the hand electrode.

- **Buzzing sound**

  You may hear a buzzing sound while operating the device.
  This is a cooling system of the device and it is not a malfunction. The buzzing sound may start and stop during use.

- **Recommended environment (Initial plate temperature of 20-30ºC)**

  When used at lower temperatures during winter, the HOT lamp may take longer to stop flashing and turn solid red (indicating it is ready to use).
  In some conditions, the safety function may automatically turn off the device. In such a case, turn the power back on and try again.

- **Plate temperature**

  Temperature of the plate will normally stabilize 90 sec* after the power is turned on.
  (Specifications: Warm plate: 42ºC, Cool plate: 12ºC)
  * At room temperature of 25ºC.
  * The time required for the plate temperature to stabilize depends on the ambient temperature, initial temperature of the plate, remaining battery, etc.
  * You may sometimes feel the Cool plate become less cool during a treatment.

- **Plate operation**

  The Warm and Cool plates will start operating at the same time power is initially turned on.
  The pink panel side is "Warm", silver panel side is "Cool".

Charging the Battery

This device can be operated using either the rechargeable battery or the AC adapter. When operating on battery power, be sure you have fully charged it first.

・This device uses Lithium ion battery. You may recharge the battery even if it is not fully discharged.
・Before using for the first time, or if you have not used the device for an extended period, fully charge the battery and make sure that the Charge Lamp has turned from flashing red to solid red indicating a full charge.

● Charging Time

Charge Time	Operating Time
3 Hours	40 Minutes


* Ambient temperature of 25°C
* Operating time is how long it can be used immediately after a full charge.

1 Insert the AC adapter plug into the adapter socket of the main unit.



2 Plug the AC adapter into the wall socket.

・Charge Lamp will flash red and the device will start charging.
・Charge Lamp will turn solid red when the charging is complete.
・The charging time may vary depending on the remaining battery power or the environment (ambient temperature).
・If you press and hold the POWER/MODE button while charging, the power turns on (CLEAN lamp turns solid green) and charging will stop.
NOTE: The battery will not charge while using the device.



3 When charging is complete, unplug the AC adapter from the wall socket and then remove the adapter plug from the main unit.

Treatment

How to hold the device

Hold the device as shown at right.
・During use, hold the device so that you do not press POWER/MODE button.
・Hold the device so that you DO NOT COVER THE VENTS with your hand.
・Be sure your HAND IS IN CONTACT WITH THE HAND ELECTRODES.
・Do not touch the plate which you are not applying to the skin as this will affect Ion output and the effectiveness of treatment.






Do not touch the other side of the plate.




Attaching cotton to the Warm plate

1. Lift the cotton holder and then rotate it towards the Warm plate side.



2. Place cotton on the Warm plate and press holder back in place.



3. Grasp both sides of the holder with your fingers and push down to lock in place (you will hear a click).



* Recommended cotton size is 70mm x 50mm or larger.

Treatment Modes

There are 5 treatment modes available.
Choose the desired treatment according to your skin condition.

Mode	Function	Ion	Use of Cotton	Skincare Product	Freq/ Duration
1 CLEAN	Skin cleaning by Ion cleansing	Cleanse	Yes	Thin lotion	2-3 times per week Maximum 10 min. per treatment
2 MOISTURE	Deliver skincare ingredients deep into skin* with Ion and moisturize the skin. (* to stratum corneum)	Infusion	Yes	Thin lotion	
			No	Facial Mask/ thicker lotion	
3 SMOOTH	Deliver skincare ingredients to skin surface and moisturize the skin.	Infusion	No	Gel or skin milk	
4 EYE CARE	Deliver skincare ingredients around the eye area with Ion and moisturize the skin.	Infusion	No	Gel or serum	
5 HOT&COOL	Improve skin elasticity and texture with Warm & Cool treatment. Can also use on neck and shoulders.	N/A	No	* Can use over makeup	Daily 10 min. per treatment

・When using modes 1-4, Ion is produced on both sides.
・Ion cleansing and iontophoresis are only effective when used together with a suitable skincare product.